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English

Dear Customer,
Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.

SUITABLE FOR

ACV602

- Anterolateral instability and a tendency to slide when seated.

ACV605

- A tendency to slide when seated.

ACV610

- Bedridden patients requiring control over their mobility, abdominal region and side in addition to bed restraint.
- Enables minor rotation changes in posture to the right or left.

FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.

If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

ACV602

- Adaptation:

Place the jacket around the thorax and the abdomen whilst the patient remains seated, then pass the waist and shoulder straps over the back of the chair, securing the patient and the back support.

Adjust the length with the sliding buckle, keeping the fasteners tightened, then pass the perineum strap between the legs, introducing them through the gap between the back support and the seat using the fastener to join the end. Finally, adjust the different straps (abdominal, perineum and shoulders) until the required securement is achieved.

The perineum strap can be used in seats in which the back support is separated from the seat, in which case the strap can be introduced between the two. Should the back support and the seat be joined together in a single piece, the perineum strap should fasten round the seat, sliding the strap underneath towards the back support where the corresponding buckle for securement and fastening is found.

ACV605

- Adaptation:

Place the belt around the patient's abdomen, fastening it to the back support using the securement straps. Pass the lower strap through the perineum zone between the thighs. Join both ends of the abdominal straps at the rear of the back support using the fastening system, pulling both straps tight until the required securement is achieved, then attach the perineum strap to the T-shape securement piece and adjust.

The perineum strap should go round the patient and slide under or through the seat towards the back support, where the corresponding buckle can be found for securement and fastening.

ACV610

- Positioning on the bed:

Extend the harness over the bed at the height of the abdominal zone with the hand clasps towards the feet.

The ends of the lower part have attachment straps with holes reinforced with metal rings, to be passed through the bed frame.

The upper part has a series of slides and straps, also with holes reinforced with metal rings, which allow us to adjust the degree of rotation permitted for the patient.

Join the frame straps and the abdominal corset strap using the magnetic fastener, introducing the metal rod and attaching with the set button. Once the abdominal corset has been adapted, use the metal rings to select the required range of rotation movement for the patient.

Prolongation straps are included for use with a wider bed.

- Positioning on the patient:

Completely extend the upper belt and lie the patient down, ensuring the strap coincides with the patient's abdominal zone. Secure the patient by overlapping one end on the other.

Introduce the rod through one of the orifices (metal rings) of the lower tab, which is protected with a strap of the same fabric, therefore preventing the loss of the rod in opening and closing movements. Then extend the upper tab, girdling it to the patient and adjusting, and introduce the rod in the required orifice. Finally, place the set button on the rod, ensuring it is perfectly in place.

Ensure all the components are properly in place and the mobility range required for the patient is correct. This system allows the patient to move his or her arms and legs, sit down and make small sideways movements. Should the immobilisation of arms and legs be required, use the components designed for this purpose.

PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.


To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the (Ref.: ACV610) symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.

RECOMMENDED CARE INSTRUCTIONS

Hand wash from time to time using warm water, pH neutral soap, and allow to dry naturally. Store in a dry place at room temperature.

KEEP OUT OF THE REACH OF CHILDREN.

This article is defined as a class I medical device.

The manufacturing date is included in the batch number which can be found on the packaging label as in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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